# User's Instructions

AH1 - General Purpose Gripping Aid

AH2 - Outdoor Gripping Aid

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Active Hands General Purpose Gripping Aids (AH1) are padded to reduce chafing. Outdoor Gripping Aids (AH2) are designed to be worn over gloves, protecting your hands from the weather and from chafing.

Your Gripping Aid has been sent to you already fastened in the gripping position, which is how the aid will look when you are wearing it.

To put on the aid please follow the following instructions whilst referring to the corresponding diagrams.

- Undo the Velcro fixing on the wrist strap.
- Undo the Velcro fixing on the gripping strap and pull Ring B through Loop B. Your aid should now appear as in Diagram 1.
- Loosen the wrist strap, but do not pull the restraint through Loop A.
- Slip the aid over your hand.
- With the help of Ring A, pull the wrist strap up tightly and fasten with the Velcro, with Loop A lying on the opposite side of the wrist to the thumb. Loop B should then be on the palm side of your hand.
- Place the object you wish to hold in your palm between the fingers and thumb, ensuring that the thumb is not bent in an uncomfortable position.
- Insert Ring B into Loop B and fold back, pulling the hand into a fist shape.
- Fasten the Velcro over the top of the fist to an appropriate level of tightness.



## Precautions

If you are unsure whether this gripping aid is suitable for you, ask your Doctor, Physiotherapist or Occupational Therapist before use.

Always follow the instructions on this leaflet when putting on and using the gripping aid.

- Active Hands gripping aids have been designed to assist with the gripping of objects which you may have been unable to use for some time. Therefore, care should be taken not to overstrain muscles (particularly in the wrists) when using heavy or awkward objects.
- Ensure the gripping aid is fastened tightly enough to hold the object firmly, so that it cannot slip from your grasp.
- Do not use the gripping aid for prolonged periods, particularly if you suffer from problems with circulation.
- Do not use the gripping aid for holding sharp objects if you suffer from blackouts.
- If the gripping aid becomes damaged, do not use.
- Do not use the gripping aid to hold hot objects or near open flames.
- When using the gripping aid to perform a new task it is advisable to have another person present.
- When using the gripping aid, remember that once it is fastened you cannot "let go" automatically.
- If holding heavy or bulky objects (e.g. weights) in the air make sure you have enough strength and control to lower them safely and carefully. Be particularly careful when lifting objects over the head.
- When undertaking exercises such as weight lifting always remember to save enough energy to undo the gripping aids.

*The use of the gripping aids may enable you to take part in pursuits which entail some measure of risk (e.g. using sharp or heavy tools, skiing, etc). Such activities should be undertaken only with the necessary supervision and recognising the inherent hazard involved.*

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