



DESCRIPTION:

A unitary device which allows a disabled person to pull a button through a button hole and pull a zipper with zero grip or dexterity (passive grip, natural tendonesis) by means of using a thumb hole.

USAGE:

Place device over the thumb, past the knuckle with the flat side toward your body. (Common for both of the following procedures)

For the zipper:

Maneuver the "Zipper tool" through the opening in the zipper tab at the bottom (closest to the zipper teeth). Place hand flat against the waist and move hand in the direction parallel to the zipper maintaining at least half of the thumb hole above the knuckle by holding the hand flat against the waist.

For the button:

This next step may be easier using the left hand, then changing to the right hand for the rest of the task. Insert the "Button tool" through the button hole from the outside of the pants.

Move the large side of the tool over the button. Rotate the opening of the tool against the top of the button until one side drops inside of the tool. Pressing against the waist, move the tool toward the center of the button while rotating until the complete button drops inside the tool.

Pull the button through the button hole by rotating while pulling.

Disengage the tool by sliding the large side of the opening toward the button, rotate and pull away from the button.



