

My Very Important Quick Start Guide

Hello, my name is Reminder Rosie. Thank you for bringing me into your home. I am thrilled to be your *personal reminder assistant* for many years to come with no conditions! Once reminders or messages are pre-recorded I will announce them in the comforting voice of a loved one to help you remember important tasks. You will hear messages over 100 feet (30m) away, accomplish the task and turn off the reminder by saying **“Reminder off”** or by pushing me down firmly. **That’s it!** Reminders can be pre-recorded in any language & announced at a specific time every day, weekly, by date or annually. (I only understand English commands).

I can help organize your day if you want me to. Just wake me up by either saying **“Hello Reminder Rosie”** or press me down firmly, and then ask: **“What day is it?”**, **“Today’s reminders”**, **“What time is it?”**, **“Help”** & more!

My brain is programmed with the latest cutting-edge voice recognition technology. I will help you coordinate all the things you need to do every day to stay healthy. To get started, invite a family member, neighbor, caregiver or health care provider to record messages. People love to hear familiar voices especially if they live alone. Please read this page & the user manual. Receive email notice of a new version & a **free 6-month extended warranty** by emailing rosie@lat.care : just type ‘Extend My Warranty’ in the subject line. Talk Soon!

Helpful Hints

- I understand normal conversational English, spoken clearly facing me and my ear (**MIC**rophone) at arms length.
- Ensure no background noise is present when setting Rosie (i.e. TV, people talking, fan, loud outside noise).
- Wake me up (trigger) **always** before saying commands; say **“Hello Reminder Rosie”** or *push* down **very firmly**.
- I cannot understand commands when I am talking. Wait about half a second after I stop talking before you do.
- I will stop talking (sleep) by pushing me down firmly anytime, once or twice. (Note two buttons on the bottom).
- I have a personality just like you & may take some practice to set reminders. Initially, you *may* need to repeat the time or date just as you would meeting a person for the first time. I may understand better when you emphasize the first syllable or first word of a 2-3 word phrase, but always in a smooth conversational voice.
- I automatically play a ‘Get Started’ introduction when plugged in without battery back-up. You can place 3 AAA quality batteries after you plug me in. Place me in an area that the user spends the most time.
- The commands I understand are printed on my back. Place the tent card beside me for quick reference.
- Ensure the user is shown how to turn off reminders; just say **“Reminder off”** during the 5-second silence between reminder playbacks or **by pressing me down very firmly**. If the user misses a reminder(s) (not turned off within 30 minutes), the REMINDER light blinks. Just wake me up to hear the missed reminder(s).
- If I malfunction, do a reset: Unplug. Remove batteries. Plug in. Set time & date. Replace batteries.(reminders stay)

Set Time and Date

- Wake me up & say **“Set time”** & follow my prompts. (Note: “Hello Reminder Rosie” phrase printed on my face).
- Wake me up again & say **“Set date”** & follow my prompts. Setting the date manually (one-time only setting with battery back-up) is recommended for *some voices*. See ‘*Setting Date Manually*’ section in the user manual.

Record, Review and Erase Reminders

- The lead caregiver pencils a Reminder List (see sample in the User Manual), the users’ important tasks with respective time/day/dates. Record up to 25 messages, 6 seconds long each. Obtain all the important reminders from family members, caregivers & healthcare providers (i.e. OT, Pharmacist, etc.). Please include both health and non-health messages & tasks, i.e. meds (how/when to take them), treatment regimens, exercise, doctor appointments, post-op care, etc. Personal reminders are just as important i.e. social, household, favorite TV shows (with channels), birthdays, and my favorite “Good night mom/dad/grandma/grandpa, I love you!”
- Ask other loved one’s (brother, sister, grandkids) to be present to record reminders. Comforting & loving voices heard everyday has shown to help people adhere to meds & other tasks effectively and over extended periods of time. I help engage the users’ ‘community of care’ by increasing everyone’s awareness of all the daily tasks.
- Now, wake-up Rosie and say **“Record reminder”**, & follow my prompts. Record reminders very loudly for louder playbacks. Remember, I take a little time getting used to. If I heard the wrong date the first time, continue one more prompt & you can tell me the full date again.
- Review or **“Erase”** saved reminders by asking me **“Play all reminders”**, or use manual mode for faster review.

Also, important: If I start talking inadvertently, say **“Cancel”** or press me down. If this is an ongoing problem, i.e. a noisy environment, turn off the TRIGGER (I won’t wake-up verbally, but I will do everything else). I will continue to announce reminders and understand commands during limited power outages even when my display is off, but I trigger only manually. I remember reminders even without batteries. You cannot change a reminder time once set; simply erase & re-record. PRINTED IN CHINA **v6**