

CUPHOLDER

Mounting Instructions

THE WRIGHT STUFF



1. Cupholder as shipped from Manufacturer.



2. Release velcro; feed strap around armrest and through the metal ring, pulling tight.



3. Position top of cupholder at the top of the armrest. When in desired position lock velcro.



4. When drink is placed in cupholder pull the cinch tight to hold drink in place.



5. Enjoy a refreshing drink with your new cupholder!



For Vertical Cupholder

Directions remain the same as above, though placement of velcro strap is designed for vertical mounting positions as can be seen in the picture to the left.

