

Made in U.S.A.

GREAT STUFF!
Scooters
Wheelchairs
For & Walkers



Diestco Manufacturing Corporation • P.O. Box 6504 • Chico, CA 95927 • 1-800-795-2392

Walker Holder for Scooters

Part No. B6214

Now with 300 yard visible reflective safety tape!!



Attaches to your scooter seatback to carry your basic style walker wherever you go!

Products Proudly Made in the U.S.A.

 **DIESTCO Mfg. Corp. • 800-795-2392 • www.diestco.com**

Walker Holder For Scooters

Attaches to the back of your scooter to carry your basic style walker. Made from high quality, water resistant polyester.

Mounting Instructions: Slip the horizontal webbing (fed through the two vertical loops) over the top of your seatback so that the bottom pocket hangs from the back of the seat. To tighten the harness, the 1" webbing utilizes a D-Ring and fastens over itself, securing with hook and loop.

After Walker Holder is secured to the seatback place your *folded* walker down into the bottom pocket. Pull the fastening straps around your walker and buckle, tighten if necessary. If your walker does not fit into the bottom pocket, either turn the walker upside down, or place only 2 or 3 of the legs into the pocket.

OTHER EQUIPMENT HOLDERS & SEATBACK BAGS AVAILABLE FOR YOUR SCOOTER OR POWERCHAIR:

Crutch Holder: Fits forearm and underarm style crutches. Goes behind your scooter or powerchair seat so you can travel without interferences.

O2 Holder: A very safe way of carrying your oxygen cylinder on the back of your scooter or powerchair.

Cane Holder: Attaches to the back of your scooter or powerchair to carry your Cane.

Deluxe Seatback Bag: Our most popular bag! Several storage compartments to separate your belongings. Wide enough to carry plenty of school and office supplies.

Monster Bag: As big as they come! Great for shopping, laundry, groceries, etc. Frame allows bag to remain open.