


AquaSense Adjustable Toilet Safety Rails, to Floor

HCPCS:



The AquaSense® Adjustable Toilet Safety Rails make it safer and easier to get on and off the toilet, it is aesthetically nice and extremely easy to assemble.

- Both lightweight and durable the rust-resistant aluminum frame supports up to 250 lbs
- Each leg includes a non-skid, non-slip tip offering security on wet surfaces
- Find added comfort and support with plastic armrests
- Simply and quickly adjust the height for proper fit and comfort. Height adjustments on this model range from 28" to 32"

PRODUCTS	Description
770-665	AqS Toilet Rail, On-Floor, 1EA,4CV,Ret
SPECIFICATIONS	Description
Warranty	Limited Lifetime
Weight Capacity	250 lb