

ELK

The ELK is a lifting cushion designed to lift from the floor.



Lifts in confined spaces



Technical Specifications

Max user weight	450 kg / 70 stone
Max height	56 cm / 22 "
Min height	Flat
Seat depth	50 cm / 20 "
Width	57 cm / 22 ½"
Total Product weight	3.6 kg / 8 lbs
Airflo Compressor	Plus

Product code: MPCA070200

Easy to use

- Lifting cushion is inflated using an Airflo Compressor with a simple hand control
- Simple to use, requiring very little training
- Suitable for use inside or out
- Can be used by one carer

Safe lifting

- Lifts the frailest person up to 70 stone / 450kg from the floor
- Exceptionally rigid and stable
- Ideal for carers of individuals with neurological disorders who fall regularly

Simple to maintain

- Clean using a proprietary non-abrasive mild liquid cleaner.
- Easy to deflate, fold away and store
- Does not need a LOLER certificate

Improves wellbeing

- Provides a safe, dignified lift for a fallen person
- Reduces the risk of musculoskeletal injury for the carer
- Can be used in confined spaces where traditional lifting equipment may not fit

