

# Handy Pillowlift

The Handy Pillowlift assists people who need help to rise from lying to sitting.



Can be used in a double bed without disturbing the non-user



## Technical Specifications

Max user weight	135 kg / 21 stone
Min height	Flat
Total length	80 cm / 31 1/2"
Width	61 cm / 24"
Total Product weight	2 kg / 4 1/2 lbs
Back/headrest length	73 cm / 28 3/4"
Backrest angle adjustment	Flat to 55 degrees
Airflo Compressor	Mk3

Product code: MPCA120200

### Easy to fit

- Sits on top of the mattress
- With no fixtures or fittings, it can be used on any size bed
- One airline connects the Airflo to the Pillowlift. Press one button to inflate

### Safe to use

- An ultra-lightweight inflating backrest raises a person from lying to sitting
- Backrest moulds around a person to provide support
- Inflation and deflation can be stopped at any time, allowing freedom to choose the most comfortable position.
- Anti-slip pillow holder

### Simple to maintain

- For a routine clean, sponge clean with a detergent solution and dry thoroughly
- The covers can be machine washed
- Spare covers can be purchased

### Improves wellbeing

- Stable, supportive and comfortable, a person can choose to be lifted on their back or side
- The Handy Pillowlift sits on top of the mattress and causes no disturbance to a second person in a double bed
- Minimises moving and handling risks so reducing the risk of musculoskeletal injury for the carer
- Its portable and lightweight which means it can be rolled up for travelling

